



BIOTEKE

Focus On Respiratory Health



FAST

3 Minutes Test
15 Minutes Results



MULTIPLE

Only 3 Drops of
Samples, Detect
multipathogen
Antigens at Once




SAFE

Only Human
AnNasal Swab
Specimens Required



SIMPLE

Easy-to-use,
No Special
Equipments Needed



WWW.BIOTEKE.CN

IVD CE3018

One-step, Detect 5 Respiratory Diseases

SARS-CoV-2

Flu A

Flu B

RSV

ADV



A


Preventive testing, caring for your health

B

Screening for suspected infection.
Auxiliary diagnosis, appropriate treatment.

C

Prognosis monitoring and confirmation of
treatment effect.



INTENDED USE & SCENARIOS



FOR SELF-TESTING

Intended Use

The kit is only used for in vitro qualitative detection of SARS-CoV-2 / Flu A / Flu B RSV / ADV antigen from human anterior nasal swab specimens.

Scenarios

- a.Home self-testing
- b.Clinics&Pharmacy
- c.Hospitals&Clinical
- d.Primary health care sites
-



ADVANTAGES



- Quick diagnostic testing, results in 15 minutes.
- Accurate results, accuracy 99%
- Easy-to-use, only nasal swab specimens needed.
- Compliance certification, Comply with EU IVDR regulations.

01 PERFORMANCE

Sensitivity > 90%
Specificity > 99%
Accuracy > 99%

Reliable

02 FAST


The pathogen can be
detected at the initial
stage of infection

Speed

03 SIMPLE

No professional
technicians and
equipment are required.

Use



PACKING SPECIFICATION









Name	Components	Loading Quantity		
Multiple Respiratory Multipathogen Antigen Test Kit		1 test/kit	5 tests/kit	20 tests/kit
	Test card	1pc	5pcs	20pcs
	Tube cap	1pc	5pcs	20pcs
	Sample extraction tube	1pc	5pcs	20pcs
	Swab	1pc	5pcs	20pcs




User manual, waste bag, QC certificates are included in the kits

TEST PROCEDURE







RESULTS INTERPRETATION



POSITIVE



NEGATIVE



INVALID